


14030, 14048, 14054, 14060 Pressure Prevention Product Drawings



14030

3 Hose Pump



- 1. Mattress Only 3 Hose- 14030MN
- 2. Mattress Cover- 14030C
Mattress Base- 14030-BASE
- 3. Bladder
w/o Laser Holes- 14030CELL
with Laser Holes- 14030-LASERCELL
- 4. "L" Connector- 14001-L
- 5. "T" Connector- 14001-T
- 6. One Way Valve- 14028V
- 7. CPR Valve- 14013
- 8. 7" Tubing- 14028T

- 9. Pump 3 Port- 14030XP
- 10. Pump Socket- 14030XP-S
- 11. Hanging Hook- 14030HK
- 12. EZ Lock Connector- 14007-3PORT
- 13. "O" Ring (3 ea)- 14007OR-3
- 14. Intake Filter- 14029FN
- 15. Power Cord- 15005PC [1]
140296G-PC [2]
Power Receptacle- 14030PR [1]
14030PR-6G [2]

[1] Three prong s/n (21)6G201210796 & before.


[2] Two Prong s/n (21)6G2010797 & later.

07/24/2023

14030, 14048, 14054, 14060 Pressure Prevention Product Drawings



14048 3 Hose Pump



- 1. Mattress Only 3 Hose- 14048MN
- 2. Mattress Cover- 14048C
- 2. Mattress Base 14048-BASE
- 3. Bladder
- 3. w/o Laser Holes- 14048-CELL
- 3. with Laser Holes- 14048-LASERCELL
- 4. "L" Connector- 14001-L
- 5. "T" Connector- 14001-T
- 6. 1 Way Valve- 14028V
- 7. CPR Plug- 14013
- 8. 7" Tubing- 14028T
- 9. Pump 3 Port- 14030XP [1]

- 10. Pump Socket- 14030XP-S
- 11. Hanging Hook-14030HK
- 12. EZ Lock Connector- 14007-3Port
- 13. "O" Ring (3 ea)-14007OR-3
- 14. Intake Filter- 14029FN
- 15. Power Cord 3 Prong Plug- 15005PC
- 15. Power Receptacle 3 Prong Plug- 14030PR
- 15. Power Cord 2 Prong Plug- 140296G-PC
- 15. Power Receptacle 2 Prong Plug-14030PR-6G

[1] Replacement pump may be 2 prong plug.


10/12/2023

14030, 14048, 14054, 14060 Pressure Prevention Product Drawings



14054

3 Hose



- 1. Mattress Only 3 Hose- 14054MN
- 2. Mattress Cover- 14054C
- 2. Mattress Base 14054BASE
- 3. Bladder
- 3. w/o Laser Holes- 14054-CELL
- 3. with Laser Holes- 14054-LASERCELL
- 4. "L" Connector- 14001-L
- 5. "T" Connector- 14001-T
- 6. 1 Way Valve- 14028V
- 7. CPR Plug- 14013
- 8. 7" Tubing- 14028T
- 9. Pump 3 Port- 14030XP [1]

- 10. Pump Socket- 14030XP-S
- 11. Hanging Hook-14030HK
- 12. EZ Lock Connector- 14007-3Port
- 13. "O" Ring (3 ea)-14007OR-3
- 14. Intake Filter- 14029FN
- 15. Power Cord 3 Prong Plug- 15005PC
- 15. Power Receptacle 3 Prong Plug- 14030PR
- 15. Power Cord 2 Prong Plug- 140296G-PC
- 15. Power Receptacle 2 Prong Plug-14030PR-6G


[1] Replacement pump may be 2 prong plug.

10/12/2023

14030, 14048, 14054, 14060 Pressure Prevention Product Drawings



14060 3 Hose



- 1. Mattress Only 3 Hose- 14060MN
- 2. Mattress Cover- 14060C
- 2. Mattress Base 14060-BASE
- 3. Bladder
- 3. w/o Laser Holes- 14060-CELL
- 3. with Laser Holes- 14060-LASERCELL
- 4. "L" Connector- 14001-L
- 5. "T" Connector- 14001-T
- 6. 1 Way Valve- 14028V
- 7. CPR Plug- 14013
- 8. 7" Tubing- 14028T
- 9. Pump 3 Port- 14030XP [1]

- 10. Pump Socket- 14030XP-S
- 11. Hanging Hook-14030HK
- 12. EZ Lock Connector- 14007-3Port
- 13."O" Ring (3 ea)-14007OR-3
- 14. Intake Filter- 14029FN
- 15. Power Cord 3 Prong Plug- 15005PC
- 15. Power Receptacle 3 Prong Plug- 14030PR
- 15. Power Cord 2 Prong Plug- 140296G-PC
- 15. Power Receptacle 2 Prong Plug-14030PR-6G

[1] Replacement pump may be 2 prong plug.

10/12/2023